



Item

AX-2000A Vacuum Mixer

Article Number

MP01001

Certification

CE Approved ISO9001:2000 Approved

Purpose

The Machine can be used to mix plasters, investments and silicones.

Description

The AX-2000A Vacuum Mixer is an earlier model, it has all the basic functions of a dental vacuum mixer and a built-in vibrator, it is one of the best selling models for its durability and lower price. It provides Hands-free and bubble-free mixing thanks to the built-in vacuum pump. This unit also features high vacuum and mixing power, and adjustable vibration.

Technical Data

Power Supply: 220V/110V 50/60Hz

Power Consumption: 150W

Motor Rotational Speed: 320rpm

Mixing Beaker: 750ml

Vacuum: 0.08~0.09MPa

Built-in Vibrator: Variable Speed

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy1@hotmail.com



Aixin Medical Equipment Co.,Ltd.China
中國天津。愛鑫醫療設備有限公司

Dimensions: 53cm x 30cm x 33cm (H W D)

Shipping Dimensions: 62cm x 38cm x 40cm (H W D)

Weight: 28kg

Parts Included

Machine×1

Manual×1

Power cord×1

Fuse×2

Mixing Beaker×1

Motor Oil×1

Filter Tap×1

Rubber Gasket ×2

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com