



**Item**

AX-4-10B Burnout Furnace (Intelligent)

**Article Number**

MW01004

**Certification**

ISO9001:2000 Approved

**Purpose**

This unit is used for the wax elimination of casting rings and preheating the crucible in dental laboratory.

**Description**

The AX-4-10B Burnout Furnace was developed in 2008. Comparing to conventional materials, the ceramic fiber furnace walls allows quick heat rise and cooling, only 15 minutes is needed for the temperature to reach 1,200°C. The furnace door is made of light-weight refractory materials. It also features automatic control system, smooth door movement and excellent thermal insulation performance. The automatic control system of the fan ensures a clean working environment.

**Technical Data**

Power Supply: 220V/110V 50/60Hz

Power Consumption: 4,000W

Temperature Rating: 1000°C

**Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com

Heating Elements: Ceramic Fiber

Capacity: Large Casting Rings x 7 or Large Casting Rings x 4 and Small Casting Rings x 4

Chamber Size: 12 × 30 × 20cm (H W D)

Dimensions: 70cm x 64cm x 47cm (H W D)

Shipping Dimensions: 74cm x 72cm x 55cm (H W D)

Weight: 65kg

---

### **Parts Included**

Machine×1

Manual×1

Thermocouple×1

Thermocouple Wire×1

### **Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com