



**Item**

AX-B3 Fine Blasting Unit

**Article Number**

DW09001

**Certification**

CE Approved ISO9001:2000 Approved

**Purpose**

This unit is used for the gentle finishing of metal frameworks, and for gross-removal of investment material and oxidation layer from clasps, crown and bridge, and rings etc.

**Description**

AX-B3 Fine Blasting Unit features large viewing window, two pencils for media of different sizes and hardness, transparent sand bottle for observing the state of sand during sandblasting. It also features strong sandblasting power, smooth and accurate operations, and low maintenance.

**Technical Data**

Power Supply: 220V/110V 50/60Hz

Power Consumption: <8W

Working Pressure: 4Kg/cm

**Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com

Handpiece: Twin Pencil

Nozzle: Replaceable

Nozzle Materials: Carbide/Ceramic

Nozzle Diameter: 1.2mm/0.8mm

Particles Size: 120µm/50µm

Optional Accessories:Dust Extractor (Reference AX-SUPER800,AX-MX800)

Dimensions: 28cm × 42cm × 40cm (H W D)

Shipping Dimensions: 36cm × 47cm × 49cm (H W D)

Weight: 18kg

---

### **Parts Included**

Manual × 1

Nozzle × 3

Lighting Tube × 1

Fuse × 2

Air Pipe × 1

Dust Pipe × 1

Aluminum oxide × 1

### **Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com