



**Item**

AX-D1 Electrolytic Polisher

**Article Number**

RW02001

**Certification**

ISO9001:2000 Approved

**Purpose**

The unit can be used for polishing cr-co and stainless steel workpieces.

**Description**

This AX-D1 Electrolytic Polisher features two water baths, high efficiency, easy operations, low maintenance, and compact design which is suited for dental laboratories.

**Technical Data**

Power Supply: 220V/110V 50/60Hz

Power Consumption: 150W

Timer:0~15min

Working Current:10A

Dimensions: 23cm x 30cm x 26cm (H W D)

**Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com



Aixin Medical Equipment Co.,Ltd.China  
中國天津。愛鑫醫療設備有限公司

---

Shipping Dimensions: 30cm x 37cm x 33cm (H W D)

Weight: 12kg

---

#### **Parts Included**

Machine × 1

Manual × 1

Power cord × 1

Fuse × 4

Alligator Clip × 1

#### **Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com