



Item

AX-J2 High Speed Alloy Grinder

Article Number

RW03003

Certification

ISO9001:2000 Approved

Purpose

The unit is for polishing and grinding all kinds of metals in dental laboratories.

Description

AX-J2 High Speed Alloy Grinder features low noise, high rotational speed, easy operations. It is equipped with safety eye shield and lighting system. The high power motor make little noise, and the high speed spindle can reach 14,000 RPM.

Technical Data

Power Supply: 220V/110V 50/60Hz

Power Consumption: 370W

Rotational speed: 14,000rpm

Optional Accessories: Dust Extractor (Reference AX-SUPER800,AX-MX800)

Dimensions: 30cm × 35cm × 29cm (H W D)

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com



Aixin Medical Equipment Co.,Ltd.China
中國天津。愛鑫醫療設備有限公司

Shipping Dimensions: 37cm × 42 cm × 36cm (H W D)

Weight: 16kg

Parts Included

Machine

Manual × 1

Power Cord × 1

Fuse × 2

Leather Belt × 1

Hexagon Wrench × 1

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com