



**Item**

AX-J3 Dental Laboratory Polishing Lathe

**Article Number**

RW03003

**Certification**

ISO9001:2000 Approved

**Purpose**

The unit is used for polishing castings in dental laboratories, e.g., alloys (Co-Cr, Ni-Cr) stainless steel, flexible dentures, bridges and porcelain teeth.

**Description**

AX-J3 Dental Laboratory Polishing Lathe features low noise, high rotational speed, easy operations, long working hours, and it is equipped with a safety dust collector hood which can be connected to a dust collector for dust-free operations.

**Technical Data**

Power Supply: 220V/110V 50/60Hz

Power Consumption: 370W

Rotational speed: 3000rpm

Optional Accessories: Dust Extractor (Reference AX-SUPER800,AX-MX800)

Dimensions: 23cm × 57cm × 17cm (H W D)

Shipping Dimensions: 32cm × 68cm × 29cm (H W D)

Weight: 13kg

---

**Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy1@hotmail.com



Aixin Medical Equipment Co.,Ltd. China  
中國天津。愛鑫醫療設備有限公司

---

#### **Parts Included**

Machine

Manual × 1

Power Cord × 1

Suction Hood × 2

Polishing Canvas Wheel × 2

#### **Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com