



Item

AX-J4 Dental Laboratory Polishing Compact Unit

Article Number

RW03004

Certification

ISO9001:2000 Approved

Purpose

The unit is used for polishing castings in dental laboratories, e.g., alloys (Co-Cr, Ni-Cr) stainless steel, flexible dentures, bridges and porcelain teeth.

Description

AX-J4 Dental Laboratory Polishing Compact Unit features low noise, high rotational speed, easy operations, and it is equipped with a suction box with light which can be connected to a dust collector for dust-free operations. The lathe, lighting system, safety eye shield and suction system make it convenient to use.

Technical Data

Power Supply: 220V/110V 50/60Hz

Power Consumption: 370W

Rotational speed: 3000rpm

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com

Optional Accessories: Dust Extractor (Reference AX-SUPER800,AX-MX800)

Dimensions: 46cm × 60cm × 50cm (H W D)

Shipping Dimensions: 53cm × 68 cm × 50cm (H W D)

Weight: 36kg

Parts Included

Machine

Manual × 1

Power Cord × 1

Fuse × 2

Polishing Canvas Wheel × 2

Dust Pipe × 1

Lighting Tube × 1

Thin Spindle × 2

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com