



Item

AX-SCA Steam Cleaner

Article Number

PF05001

Certification

CE Approved ISO9001:2000 Approved

Purpose

This unit is for cleaning castings and instruments, removal of grease, sticky wax, pickling agents, contaminants and adhesives.

Description

AX-SCA Steam Cleaner has both steam-gun-controlled and button-controlled steaming methods, it features low water, over-pressure, over-heat protection and forced heat function. It also has high heating power, quick steam generating power, high efficiency, strong and stable cleaning power.

Technical Data

Power Supply: 220V/110V 50/60Hz

Power Consumption: 1600W

Max. Pressure: 6Bar (default)

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy1@hotmail.com



Aixin Medical Equipment Co.,Ltd. China
中國天津。愛鑫醫療設備有限公司

Working Pressure: 4~5Bar

Water Tank Capacity:5.5L

Water Injections: Manual

Dimensions: 54cm × 38cm × 38cm (H W D)

Shipping Dimensions: 59cm × 45cm × 44cm (H W D)

Weight: 24kg

Parts Included

Machine

Manual × 1

Power Cord × 1

Self-lock Buttonl × 1

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy@hotmai.com