



Item

AX-SCB Steam Cleaner

Article Number

PF05001

Certification

CE Approved ISO9001:2000 Approved

Purpose

This unit is for cleaning castings and instruments, removal of grease, sticky wax, pickling agents, contaminants and adhesives.

Description

AX-SCB Steam Cleaner is the updated version of the AX-SCA Steam Cleaner. Users can use tap water instead of distilled water which reduces production costs. The boiler, casing, buttons, indicator lights, and internal parts are made of stainless steel which significantly increases its lifespan and reliability. This unit has both steam-gun-controlled and button-controlled steaming methods, it features low water, over-pressure, over-heat protection and forced heat function. It also has high heating power, quick steam generating power, high efficiency, strong and stable cleaning power.

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com

Technical Data

Power Supply: 220V/110V 50/60Hz

Power Consumption: 1600W

Max. Pressure: 6Bar (default)

Working Pressure: 4~5Bar

Water tank capacity: 7.5L

Water Injections: Manual

Dimensions: 57cm × 38cm × 41cm (H W D)

Shipping Dimensions: 63cm × 45cm × 49cm (H W D)

Weight: 36kg

Parts Included

Machine

Manual × 1

Power Cord × 1

Self-lock Buttonl × 1

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy@hotmai.com