



**Item**

AX-YD Denture Injection System

**Article Number**

RW05001

**Certification**

ISO9001:2000 Approved

**Purpose**

This unit is for making flexible removable partial dentures, it includes a separate manual press and furnace.

**Description**

AX-YD Denture Injection System features easy operation, long lifespan, low maintenance, and high success rate.

**Technical Data**

Power Supply: 220V/110V 50/60Hz

Power Consumption: 400W

Heating Temperature: 287±10°C

Thermocouple: K Type

Dimensions: 77cm x 55cm x 55cm (H W D)

Shipping Dimensions: 32cm x 65cm x 35cm (H W D)

**Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com

Weight: 45kg

---

#### **Parts Included**

Machine × 1

Manual × 1

Fuse × 2

Denture Flask × 1

Cartridge (8g) × 5

Cartridge (19g) × 5

Silicon Oil × 1

Brush × 1

Hexagon Wrench × 1

Copper Gasket × 2

Charge-in Tube × 1

Discharging × 1

#### **Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com