



#### Item

AX-Z2 Dental Vibrator

#### Article Number

MP02001

#### Certification

ISO9001:2000 Approved

#### Purpose

This unit be used for reducing bubbles while pouring plasters or silicones, ensuring the flawless results of duplicating moulds.

#### Description

The AX-Z2 Dental Vibrator features heavy duty, easy operations and adjustable vibration, and the vibration can be customized.

#### Technical Data

Power Supply: 220V/110V 50/60Hz

Power Consumption: 100W

Max. Load: 2kg

Dimensions: 10cm x 18cm x 22cm (H W D)

Shipping Dimensions: 16cm x 26cm x 29cm (H W D)

Weight: 5.5kg

---

#### Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy1@hotmail.com



Aixin Medical Equipment Co.,Ltd. China  
中國天津。愛鑫醫療設備有限公司

---

### **Parts Included**

Machine

Manual×1

Power Cord×1

Fuse×2

### **Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com