



Item

EM-BX1 Dental Sandblasting Cabinet

Article Number

PF03001

Certification

ISO9001:2000 Approved

Purpose

This unit eliminates the need for cleaning up the working area for dental technicians, it features dust free operations, improving the working qualities and working environment significantly.

Description

The EM-BX1 Dust Cabinet features wide operating view, cold light system, wide hand inlet, dust collection system and wooden arm rests.

Technical Data

Power Supply:220V/110V 50/60Hz

Power Consumption:<8W

Arm Rest:Removable

Optional Accessories:Dust Extractor (Reference AX-SUPER800,AX-MX800)

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com



Aixin Medical Equipment Co.,Ltd. China
中國天津。愛鑫醫療設備有限公司

Dimensions:26cm×34cm×39cm (H W D)

Package Size:31cm×39cm×44cm (H W D)

Weight:15kg

Parts Included

Machine × 1

Manual × 1

Power cord × 1

Lighting Tube × 1

Fuse × 2

Arm Rest × 2

Dust Pipe × 1

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com