



**Item**

EM-DC1 Die Separating Unit

**Article Number**

MP06001

**Certification**

ISO9001:2000 Approved

**Purpose**

The Machine is used in dental labs, it is the ideal equipment for cutting dental models. This machine can replace the traditional handsaws.

**Description**

EM-DC1 Die Separating Unit features high speed, low noise, stable rotation, and angle-adjustable working platform. The highly reinforced tungsten saw blade, which is imported from Germany(DFS), will provide excellent cutting results.

**Technical Data**

Power Supply: 220V/110V 50/60Hz

Power Consumption: 50W

Rotational Speed: 2,800RPM

Disc Size:80mm/16mm/0.25mm (OD/ID/T)

**Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com



Aixin Medical Equipment Co.,Ltd. China  
中國天津。愛鑫醫療設備有限公司

---

Optional Accessories: Dust Extractor (Reference AX-SUPER800,AX-MX800)

Dimensions: 34cm x 37cm x 47cm (H W D)

Shipping Dimensions: 40cm x 43cm x 53cm (H W D)

Weight: 20kg

---

### **Parts Included**

Machine

Manual × 1

Power Cord × 1

Fuse × 2

Saw blade × 1

### **Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com