



Item

EM-DC2 Die Separating Unit

Article Number

MP06002

Certification

ISO9001:2000 Approved

Purpose

The Machine is used in dental labs, it is the ideal equipment for cutting dental models. This machine can replace the traditional handsaws.

Description

EM-DC2 Die Separating Unit features high rotational speed, constant RPM, and adjustable model holder. The built-in lighting system and positioning laser beam makes it even more accurate. The diamond plaster disc from DFS Germany ensures excellent product quality.

Technical Data

Power Supply: 220V/110V 50/60Hz

Power Consumption: 110W

Rotational Speed: 3,500RPM

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com



Aixin Medical Equipment Co.,Ltd.China
中國天津。愛鑫醫療設備有限公司

Disc Size:80mm/16mm/0.25mm (OD/ID/T)

Optional Accessories: Dust Extractor (Reference AX-SUPER800,AX-MX800)

Dimensions: 42cm x 30cm x 37cm (H W D)

Shipping Dimensions: 50cm x 38cm x 45cm (H W D)

Weight: 30kg

Parts Included

Machine

Manual × 1

Power Cord × 1

Saw blade × 1

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com