



Item

EM-LT1 Dental Workstation (Single)

Article Number

MP07201

Certification

CE Approved ISO9001:2000 Approved

Purpose

This unit is used in dental laboratories for removal of dust resulting from machining.

Description

The EM-LT1 frame is made of steel, and the tabletop is made of high density fireproof board. It also features florescent lighting system, stainless steel plate, arm rests, retrievable air nozzle, dust extraction system, storage drawers, etc. It also features florescent lighting system, stainless steel plate, arm rests, retrievable air nozzle, dust extraction system, and storage drawers. The adjustable lighting system enhances the lighting effect.

Technical Data

Power Supply: 220V/110V 50/60Hz

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com

Power Consumption: 800W

Air Volume: 172m³/h

Max. Pressure: -16kPa

Noise Level: <68dB

Dimensions: 85cm x 60cm x 120cm (H W D)

Shipping Dimensions: 102cm x 71cm x 127cm (H W D)

Weight: 120kg

Parts Included

Adjustable 3-D Lighting System×1

Air Nozzle×1

Arm Rest×2

Backup Power Socket×1

Dust Extraction System×1

Tools Drawer×2

Storage Drawer×3

Pan×3

Acrylic Resin Tabletop×1 (Optional)

Hettich Drawer Slide (Germany)×10 (Optional)

Gas System×1 (Optional)

Technician Stool×1 (Optional, please refer to C450)

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com