



#### **Item**

EM-LT2 Dental Workstation (Single)

#### **Article Number**

MP07202

#### **Certification**

CE Approved ISO9001:2000 Approved

#### **Purpose**

This unit is used in dental laboratories for removal of dust resulting from machining.

#### **Description**

The EM-LT2 Workstation was developed in 2007, the light arm and table frame have been modified to make it more stable. The table frame is made of steel, and the tabletop is made of high density fireproof board. It also features florescent lighting system, stainless steel plate, arm rests, retrievable air nozzle, dust extraction system, storage drawers, etc. It also features florescent lighting system, stainless steel plate, arm rests, retrievable air nozzle, dust extraction system, and storage drawers.

#### **Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com

### Technical Data

Power Supply: 220V/110V 50/60Hz

Power Consumption: 800W

Air Volume: 172m<sup>3</sup>/h

Max. Pressure: -16kPa

Noise Level: <68dB

Dimensions: 85cm x 60cm x 120cm (H W D)

Shipping Dimensions: 102cm x 71cm x 127cm (H W D)

Weight: 120kg

---

### Parts Included

Fluorescent Lighting System×2

Air Nozzle×2

Arm Rest×4

Backup Power Socket×2

Dust Extraction System×2

Tools Drawer×4

Storage Drawer×0

Pan×6

Acrylic Resin Tabletop×1 (Optional)

Hettich Drawer Slide (Germany)×8 (Optional)

Gas System×2 (Optional)

Technician Stool×1 (Optional, please refer to C450)

### Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com