



Item

Magnifier Bench Lamp

Article Number

MP08002

Certification

CE Approved ISO9001:2000 Approved

Purpose

The magnifier desk Lamp be use for dental laboratory,jewellery,electronic devices and components,Watch and Clock Repair, and so on.

Description

The Supervision Magnifier Lamp be use for all type dental workstation.You can iinstall the Magnifier Lamp on your own bench by yourself also. 5x magnification, swing arm, clamp mount, great for any detailed work, perfect for porcelain sculpting, waxing, metal, and finishing.

Technical Data

Power Supply: 220V/110V 50/60Hz

Power Consumption: W

Lens Material: glass

Lens size: 7.5 "x6.2"

Diopter: 3 diopter/5 diopter

Overseas Marketing:

Add: 6B-101,No.15 Rongyuan Road,Huayuan New-Tech Industrial Park,Tianjin,China,300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com



Aixin Medical Equipment Co.,Ltd. China
中國天津。愛鑫醫療設備有限公司

Light source: T5 28W energy saving compact fluorescent lamp

Standard mount: desktop clamp

Optional mount: table base, floor stand

Colour available: white, black

Dimensions: cm x cm x cm (H W D)

Shipping Dimensions: cm x cm x cm (H W D)

Weight: 5kg

Overseas Marketing:

Add: 6B-101, No. 15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy1@hotmail.com